

Standard Anatomy Dissection Table with Creased Top

Mopec's Creased Top Dissection Tables are thoughtfully engineered to optimize cadaver dissection while prioritizing efficiency and cleanliness. Featuring a creased marine edge, these tables prevent fluid spills and ensure a smooth, rounded perimeter for effortless cadaver placement and removal from either side.



HA200

SPECIFICATIONS

Construction

- 1.5" marine edge for ease of transfers
- Includes drain hole with PVC drain valve
- (4) heavy duty 8" casters with (3) total lock and (1) directional lock
- (4) Post sockets included
- Creased, non-removable top
- Built-in tapered slope
- Full-length lower shelf

Materials

- 304 stainless steel construction with a polished No. 4 finish

Dimensions

- 86" L x 30" W x 36.8" - 37.5" H (sloped top)

Capacity

- 350 lbs

Options & Accessories

- HO005: Post leg-holders
- HO006: Book holder
- HO030: Instrument tray holder

Related/Alternate/Obsolete Product Model Numbers

Model No.	Description	Status	Notes
600023	Standard Dissection Table	Replaced	Use HA200
HA210	Standard Dissection Table	Replaced	Use HA200

MOPEC
GROUP

Mopec
ELEVATING ANATOMY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI