

27" Portable Side-Loading Cadaver Rack

Mopec body racks are made from strong stainless steel to ensure long-lasting use and easy cleaning in labs and post-mortem settings. They fit smoothly into lab setups with features like rollers and locking trays for safe and easy movement. Built for heavy loads and daily use, Mopec racks offer reliable quality you can count on.



IW27301



IW27401

SPECIFICATIONS

Construction

- Fits 27" trays (not included)
- 10 full-width rollers per tier
- 4 tray stops on each tier with rubber bumpers
- Positive-locking tray stops on each tier

Materials

- Thick-wall 304 stainless steel construction with a 180-grit finish

Dimensions

- IW27301: 83.0 L x 27.8" W x 64.7" H
Tier heights: 23.6", 42.6", 59.6"
- IW27401: 83.0" L x 27.8" W x 72.6" H
Tier heights: 17.6", 34.6", 51.6", 68.6"

Capacity

- IW27301: 1,125 lbs
- IW27401: 1,500 lbs

Compatibility

- Fits JC027 and JC027-D 27" trays
- JD501, JD5000, 600040DS, M690E and M691ES scissor lifts
- JD450, JD451, M677, and M678 forklifts, but must be used with JD460 reversible forks and with M600H-27, M600H-SC, M601H-27, and M601H-ES27 roller pallets

Options & Accessories

- Other sizes available

Note:

- Trays sold separately

Variations

Model No.	Description	Dimensions
IW27301	Three Tier Portable Side-Loading Cadaver Rack	83.0" L x 27.8" W x 64.7" H
IW27401	Four Tier Portable Side-Loading Cadaver Rack	83.0" L x 27.8" W x 72.6" H

MOPEC

GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING



mopecgroup.com
300.362.8491
Low, Madison Heights, MI