

27" Stationary Side Loading Body Rack

Constructed with heavy-duty stainless steel, this rack allows for easy side-loading of cadaver trays, optimizing workflow and ensuring long-term structural integrity. Its stationary design provides stable, reliable storage in high-demand environments.



IS27301



IS27401 TRAY NOT INCLUDED

SPECIFICATIONS

Construction

- Stationary and side loading
- Suitable for refrigerated trailers or other cooled structures
- (10) rollers and (4) rubber tray stops per tier
- Positive-locking tray latch
- Heavy duty adjustable leveling footpads

Materials

- 11GA 304 stainless steel tube construction with a 180-grit finish

Dimensions

- IS27301: 83" L x 29" W x 58.3" H
Tier Heights: 17.25", 36.25", 53.25"
- IS27401: 83" L x 29" W x 70.3" H
Tier Heights: 15.25", 32.25", 49.25", 66.25"

Capacity

- 375 lbs per tier

Performance

- Accommodates JCo27 and JCo27-D 27" trays
- Designed for loading using a JD4000, JD5000, and M690ES, or forklift and M600H-27 roller pallet

Options & Accessories

- Other sizes available

Note:

- Tray for reference only. Trays sold separately.

Variations

Model No.	Description	Dimensions
IS27301	Three Tier Stationary Side Loading Body Rack	83" L x 29" W x 58.3" H
IS27401	Four Tier Stationary Side Loading Body Rack	83" L x 29" W x 70.3" H

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI