

Autopsy Carrier Chassis

Mopec's Autopsy Carrier Chassis provides a strong, stable foundation for autopsy procedures, supporting seamless transport and positioning within the morgue or pathology suite. Built for durability and smooth maneuverability, it offers a reliable performance and compatibility with the T3603 Body Tray. Its sturdy construction and user-focused design help enhance workflow efficiency and handling safety.



600017

SPECIFICATIONS

Construction

- (4) locking caster wheels
- For use with T3603 Plastic Body tray (not included)
- Underside of T3603 tray mates to top of the carrier to securely hold the tray in place without the use of mechanical fasteners

Materials

- 11ga 1.5" SQ 304 Stainless steel frame with a 180-grit polish finish
- Flame-retardant hard plastic carrier top

Dimensions

- 80.5" L X 26" W X 32.7" to 34.3" H (lower to higher end for drainage)

Capacity

- 500 lbs

Compatibility

- T3603 (Plastic Body Tray)

Note:

- Trays sold separately.

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI