

# Recirculating Down-Back Draft Workstation

Mopec's BF400 Workstation is a compact, self-contained countertop dual-draft workstation designed for small labs performing dissection, trimming, and staining procedures. It has a shelf on top with a light beneath it, and an adjustable speed blower system. It comes with potassium permanganate filters, with optional charcoal replacement filters. The front access panel makes the filters easy to replace. The blower system filters the air and recirculates it back into the room.



BF400

## SPECIFICATIONS

### Construction

- Includes (1) BF400 Counter down-back draft workstation (recirculating)
- Shelf on top with a light beneath and an adjustable speed blower system
- Includes potassium permanganate filters with optional charcoal replacement filters

### Materials

- 304 Stainless steel, electric blowers

### Dimensions

- 21" L x 36" W x 19" H

### Utilities / Power

- 110/120 V AC/60HZ (USA)
- 10-foot hospital-grade cord

### Performance

- Filters air and recirculates it back into the room

### Options & Accessories

- BF450 (facility-vented version of this workstation)
- BF008 Charcoal filter for BF400 Qty. of 2 Req'd
- BF007 Potassium Permanganate filter for BF400 Qty. of 2 Req'd

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI