

Counter Down-Back Draft Workstation

Mopec's BF450 Workstation is a compact, self-contained countertop dual-draft workstation designed for small labs performing dissection, trimming, and staining procedures. It has a shelf on top with a light beneath it and has both downdraft and backdraft functions. It connects to your facility ventilation system, sucking harmful fumes out of your work area and right out of the building. Engineered with all-stainless-steel construction, the BF450 features Mopec's Dual-draft ventilation system for safe and effective fume removal.



BF450

SPECIFICATIONS

Construction

- Includes (1) BF450 counter down-back draft workstation (recirculating)
- Includes a shelf with light beneath

Materials

- 304 Stainless steel

Dimensions

- 21" L x 36" W x 19" H

Utilities / Power

- 110/120 V AC/60HZ (USA)
- 10-foot hospital-grade cord
- 18" x 4" rectangular vent duct connection

Options & Accessories

- BF400 (Recirculating version of this workstation)
- BF008 Charcoal filter for BF400 Qty. of 2 Req'd
- BF007 Potassium Permanganate filter for BF400 Qty. of 2 Req'd

Performance

- Both downdraft and backdraft functions
- Connects to facility ventilation system
- Sucks harmful fumes out of work area and out of building

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI