

White Fume Ventilation Adapter

Mopec's Fume Ventilation Adapter is meant so you can connect your fume extractor to the building ventilation system, with the building creating suction and removing the fumes from your working area. The BF720 fume ventilation adapter is made to go on top of fume extraction modules.



BF720

SPECIFICATIONS

Construction

- Made to go on top of the BF724 through BF760 fume extraction modules
- Connectable to building ventilation system
- Includes (1) BF720 ventilation adapter

Materials

- Powder Coated Steel

Dimensions

- Squared bottom: 12.2" L x 12.2" W x 7" H
- 5.990" Diameter connection for a 6" flexible duct

Compatibility

- Fits BF724, BF736, BF748, BF760, fume extraction modules

Options & Accessories

- BF724, BF736, BF748, BF760, fume extraction modules
- PM0603 vent duct (sold separately by the foot)
- PB0378 clamp for vent duct (sold separately, individually)

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI