

31" Autopsy Cart

Mopec's autopsy carts are made of stainless steel and are a durable design to support bodies safely and efficiently during examinations. Built with a sturdy frame and heavy duty wheels, it allows staff to transport remains with ease while maintaining a sanitary and controlled work environment.



DA031

SPECIFICATIONS

Construction

- Fits GA031 31" tray (sold separately)
- (3) 8" heavy duty casters with total lock capabilities
- (1) 8" heavy duty caster with directional lock capabilities
- 3" Incline to allow for the flowing of liquids and debris

Materials

- 16ga 1.5" Square 304 stainless steel crossbars
- 11ga 1.5" Square 304 stainless steel legs
- (4) 8" heavy duty casters

Dimensions

- 77.5" L X 31" W X 34.4" H to 37.3"

Capacity

- 500 lbs

Compatibility

- GA031 31" tray with sloped sides (sold separately)

Note:

- Trays sold separately.

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI