


Aluminum False Frame and Cover

The Mopec Aluminum False Frame is a lightweight, removable frame designed to fit onto compatible Mopec cadaver carriers. Made of aluminum and paired with a vinyl cover, it gives the carrier the appearance of a standard gurney during transport. This frame helps ensure discretion and dignity, while remaining easy to install or remove as needed



DC410

SPECIFICATIONS

Construction

- The aluminum false frame can be lifted and removed from the cadaver carrier for ease of cadaver transport
- The cadaver cover provides discrete transport of the cadaver within a workplace

Materials

- Lightweight aluminum frame
- Vinyl fabric cover (JZ010, white vinyl)

Dimensions

- Frame: 81.1" L X 40.6" W X 18" H
- Cover: 81" L X 40" W X 20" D

Compatibility

- DC400 Bariatric autopsy cart
- JZ010 White vinyl cover (included)

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI