

# Hydraulic Downdraft Dissection Table

Mopec's Hydraulic Downdraft Dissection Table provides a sturdy, easy-to-use surface for dissection work. It features downdraft ventilation to help control odors and airborne contaminants, along with hydraulic height adjustment for user comfort. Built with durable materials, it's well suited for labs, schools, and forensic settings.



HB300

## SPECIFICATIONS

### Construction

- Exhaust requirements: 350–500 CFM
- (4) heavy duty 8" casters
- Raise or lower the cart using the dual sided mounted foot controls
- Enable fifth wheel steering, neutral, or full brakes with foot controls
- The vented stainless-steel table allows fumes to be pulled down and away from the worksurface

### Materials

- Vented 304 stainless steel table with No.4 polish and drainage outlet
- Hydraulic Stryker base

### Dimensions

- 86" L (89" including duct) X 30" W X 30.25" to 43.75" H

### Capacity

- 500 lbs

### Utilities

- Recommended airflow, adjustable between 350 and 500 CFM.
- Duct stub connection: 6" round
- Ball valve with 1" flexible drain

### Options and Accessories

- HO001: Leg Rest Bridge

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI