

Hydraulic Covered Cadaver Carrier with Cover

The Mopec Hydraulic Covered Cadaver Carrier with a cover is a durable, discreet transport cart featuring a stainless-steel, non-removable tray top, lightweight frame, and opaque cover. Its foot-operated hydraulic lift allows easy height adjustment, while large locking casters and a 5th-wheel design provide smooth, controlled maneuverability for safe and dignified cadaver transport.



JA600



JA614

SPECIFICATIONS

Construction

- Non-removable tray top
- Raise or lower the cart using the dual sided mounted foot controls
- Removable false frame
- Enable fifth wheel steering, neutral, or full brakes with foot controls
- (4) heavy duty 8" casters

Capacity

- 500 lbs

Materials

- 304 Stainless steel carrier top with a No. 4 polish finish (welded to cart, non-removable)
- The false frame is made of lightweight aluminum
- Hydraulic Stryker base
- Cover
 - JA600: White rubber-backed blackout fabric
 - JA614: Black elastane fabric

Dimensions

- 81" L X 30" W X 21.25"-34"H (tray top) 37.25"-50" H (with false frame on top)
- 81.7" L X 30.7" W max including bumpers

Compatibility

- JA600: JA610 (Cadaver False Frame) and JZ005 (White Cadaver Cover) Both included and also sold separately.
- JA614: JA610 and JZ029 (Black Cadaver Cover) Both included and also sold separately.

Variations

Model No.	Description	Maximum Dimensions
JA600	Hydraulic Carrier with White Cover	81.7"L x 30.7"W x 50"H
JA614	Hydraulic Carrier with Black Cover	81.7"L x 30.7"W x 50"H

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI