

Hydraulic Cadaver Carrier

The Mopec Hydraulic Covered Cadaver Carrier with white cover is a durable, discreet transport cart featuring a stainless-steel, non-removable tray top, lightweight frame, and opaque white cover. Its foot-operated hydraulic lift allows easy height adjustment, while large locking casters and a 5th-wheel design provide smooth, controlled maneuverability for safe and dignified cadaver transport.



JA620

SPECIFICATIONS

Construction

- Non-removable tray top
- Raise or lower the cart using the dual sided mounted foot controls
- Enable fifth wheel steering, neutral, or full brakes with foot controls
- (4) heavy duty 8" casters

Materials

- 304 Stainless steel carrier top with a No. 4 polish finish (welded to cart, non-removable)
- Hydraulic Stryker base

Capacity

- 500 lbs

Dimensions

- 81" L X 30" W X 21.25" to 34" H (tray top) 37.25" to 50" H (with false frame on top)
- 81.7" L X 30.7" W max length and width including bumpers

Compatibility

- JA610: Cadaver False Frame
- JZ005: White Vinyl Cadaver Cover

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI