


# Roller Pallet Assembly, Various Sizes

Mopec's Roller Pallet Assembly is engineered to provide smooth, controlled lifting and effortless transfer of heavy loads. For hydraulic forklifts, its full-roller surface ensures safe, efficient movement of trays, racks, or remains. Built with durable materials and designed for ergonomic handling.



M600H-23



M600H-T14



M600H-27

## SPECIFICATIONS

### Construction

- For use with JC Series trays and GA tray (not included)
- (10) Nylon rollers to allow for smooth movement of trays on the pallet
- Safety Latches with Bumpers on either side of the assembly to secure the tray in place during moving.

### Materials

- 11ga 1.5" SQ 304 stainless tubing construction with a 180-grit finish
- Nylon Rollers

### Dimensions

- M600H-23: Assembly: 81.8" L X 26.5" W X 5.0" H
- M600H-27: Assembly: 81.8" L X 30.5" W X 5.0" H
- M600H-T14: Assembly: 85" L X 35.3" W X 5.0" H
- Fork Pockets: 2.0" T X 4.0" W X 0.125" thick

### Capacity

- 500 lbs

### Compatibility

- M600H-23: JC023 and JC023D trays
- M600H-27: JC027 and JC027D trays
- M600H-T14: Ga031 tray
- All: M-678 Hydraulic Lift JD405, JD450, and JD451 Forklifts with JD460 Forks
- Not for use with hand crank lifts. For 23 and 27, use M600C-23 or M600C-27

### Note:

- Trays sold separately.

## Variations

Model No.	Description	Dimensions
M600H-23	Roller Pallet Assembly for 23" Trays	81.8" L X 26.5" W X 5" H
M600H-27	Roller Pallet Assembly for 27" Trays	81.8" L X 30.5" W X 5" H
M600H-T14	Roller Pallet Assembly for 31" Trays	85" L X 35.3" W X 5" H

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI