



# Roller Pallet Assembly, End/Side Loading

The Mopec Roller Pallet Assembly is a sturdy platform designed for efficient end or side loading of cadaver trays. Its heavy-duty roller system ensures smooth movement. Built with a durable steel frame, it offers reliable, long-lasting performance for mortuary, pathology, and laboratory environments.



**M601H-ES-23**



**M601H-ES-27**

## SPECIFICATIONS

### Construction

- For use with JC Series trays (not included)
- (10) multi-directional rollers
- Safety latches with bumpers on the ends of the assembly and safety side flip latches to secure the tray in place during moving.

### Materials

- 11ga 1.5" SQ 304 stainless tubing construction with a 180-grit finish

### Dimensions

- M601H-ES-23: Assembly: 80.5" L X 23.5" W X 1.5" H
- M601H-ES-27: Assembly: 80.5" L X 27.5" W X 6" H
- Fork Pockets: 2.0" T X 4.0" W X 0.125" thick

### Capacity

- 500 lbs

### Compatibility

- M601H-ES-23: JC023 and JC023D trays
- M601H-ES-27: JC027 and JC027D trays
- Both: M-678 Hydraulic Lift JD405, JD450, and JD451 Forklifts with JD460 Forks
- Note: Not for use with hand crank lifts

### Note:

- Trays sold separately.

### Variations

Model No.	Description	Dimensions
<b>M601H-ES-23</b>	Roller Pallet Assembly for 23" Trays	80.5" L X 23.5" W X 1.5" H
<b>M601H-ES-27</b>	Roller Pallet Assembly for 27" Trays	80.5" L X 27.5" W X 6" H

**MOPEC**  
GROUP



[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI