


Molded Plastic Carrier Chassis Tray

Built with durable materials and thoughtful design, Mopec's Body Tray series ensures respectful handling and secure transport of cadavers in both routine and high-demand situations. Easy to clean and built to last, they support the highest standards in hospital and forensic autopsy settings, making them a dependable choice for pathology professionals in maintaining an effective and efficient work environment.



T3603

SPECIFICATIONS

Construction

- Hand holes on either end for carrying/lifting the tray

Materials

- Hi-tech molded fire-retardant plastic, Kydex Acrylic/PVC

Dimensions

- 82" L x 31" W x 3.6" H

Capacity

- 500 lbs

Compatibility

- IJ23 racks
- JD5000 & JD5001 cadaver lifts
- M690ES-23
- JD4000
- JD420 roller pallet
- JD305, JD405, & JD710

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI