

# Covered Cadaver Scissor Lift Carrier

Mopec's Covered Cadaver Scissor Lift Carrier is a durable, easy-to-maneuver transport solution designed for safe and discreet body handling. Featuring a smooth scissor-lift mechanism, sturdy construction, and a protective cover for privacy, it allows effortless loading, unloading, and transport while maintaining sanitary conditions in medical, forensic, and mortuary settings.



600040-DS-23



600040-DS-27

## SPECIFICATIONS

### Construction

- End-Loading only
- 16" to 56.8" elevation
- Raise or lower the carrier using the handheld controller
- Removable aluminum tubing canopy and black vinyl cover
- (3) heavy duty 6" swivel locking casters
- (1) heavy duty 6" directional locking caster
- (5) galvanized stainless-steel rollers
- Emergency stop button
- (2) locking latch mechanisms at the head and foot of the carrier

### Capacity

- 500 lbs

### Materials

- 11 Gauge square 304 stainless steel carrier frame with a No. 4 polish finish
- Galvanized stainless-steel rollers
- The removable canopy frame is made of aluminum tubing
- Black vinyl removable canopy cover

### Dimensions

- **DS-23:** 81.5" L X 26.8" W X 16" to 56.8" H (72.7" H with false frame)
- **DS-27:** 81.5" L X 30.8" W X 16" to 56.8" H (72.7" H with false frame)

### Utilities / Power

- Rechargeable battery cord accepts any standard 120V/1hp/60V outlet
- 24V High-capacity battery capable of thrust up to 6000N

### Compatibility

- **DS-23:** JC023 and JC023D
- **DS-27:** JC027 and JC027D
- Note: Meant for 3 tier racks, check your application. A 3-tier rack inside a cooler may require a 5-tier lift if the lift is on the floor outside and the rack is inside the cooler

## Variations

Model No.	Description	Dimensions
<b>600040-DS-23</b>	Scissor Lift Carrier for 23" Trays	81.5"L x 26.8"W x 56.8"H
<b>600040-DS-27</b>	Scissor Lift Carrier for 27" Trays	81.5"L x 30.8"W x 56.8"H

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI