

Maestro Monitor and Keyboard Track, Left or Right Side



AX117L



AX117R

SPECIFICATIONS

Construction

- Designed to give additional space on the grossing station for mounting monitors and keyboard arms
- Attaches on the same side as the touchscreen panel of your Maestro or Maestro Encore grossing station

Dimensions

- 2.8" W x 1.6" L x 38.5" H

Capacity

- 40 lbs

Materials

- 304 stainless steel sheetmetal and hardware

Compatibility

- Fits all Mopec Maestro and Maestro Encore grossing stations
- Supports AX032, AX034 and AX070 keyboard and monitor arms
- Not compatible with AX031 pegboard wing when used on same side

Variations

Model No.	Description	Dimensions
AX117L	Monitor and Keyboard Track, Left Side	2.8" W x 1.6" L x 38.5" H
AX117R	Monitor and Keyboard Track, Right Side	2.8" W x 1.6" L x 38.5" H

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI