

Down-Back Draft Grossing Station

Mopec's Down-Back Draft Grossing Station is designed to support a safer and more comfortable environment during surgical grossing procedures. The workstation helps manage fumes and odors at the work area while maintaining an ergonomic, user-focused design. Built for durability and everyday use, it is well suited for the demands of busy pathology labs.



BF550

SPECIFICATIONS

Construction

- Facility Ventilation pulls between 300–350CFM
- 60" High pressure hand spray hose and 1.5" drainpipe
- Top mounted, watertight G.F.C.I. Outlet
- 18" Standard Fluorescent Replacement Bulb with Lighted On/Off Switch

Materials

- 18 Gauge 304 Stainless Steel with a No. 4 Finish casing
- 16 Gauge 304 Stainless Steel with a No.4 Finish frame

Dimensions

- 21" L X 36" W X 58.7" H
- Workstation is at 37.3" H
- 18" x 4" rectangular duct

Utilities / Power

- Hospital Grade Line Cord 115v/1ph/60Hz (8' long)

Options & Accessories

- BF003 Particle Filter (1 included)
- BF005 Potassium Permanganate Filter (1 included)
- BF101 Replacement Pad (1 pack of 3 included)
- BF500 for recirculating version of this table

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI