

# 31" Autopsy Cart, Chassis Only

Mopec's autopsy carts are made of stainless steel and are a durable design to support body trays during examinations. Built with a sturdy frame, options for elevation, and locking casters, it allows staff to transport remains with ease while maintaining a sanitary and controlled work environment.



DE031

## SPECIFICATIONS

### Construction

- (3) 8" heavy duty locking carrier wheels
- (1) 8" heavy duty swivel lock wheel
- (10) Nylon rollers
- 1" Incline between the head and foot of the cart, for drainage
- (2) End cart holders

### Materials

- 11ga 1.5" SQ 304 Stainless steel frame with a No.4 polish finish
- Nylon Rollers
- Rubber Tray Bumpers

### Capacity

- 500 lbs

### Dimensions

- 79.7" L X 32" W X 35.8" to 36.8" H (1" slope for drainage)

### Compatibility

- Fits GA031 trays

### Note:

- Trays sold separately

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI