


REV	DESCRIPTION	DATE	BY	APP
A	INITIAL RELEASE	00-XXX-19	XXXXX	MM
B	SEE EDN	1-7-26		DA




(4) REMOVABLE
GRID PLATES

MOPEC LID ASSIST
HINGED COVER W/ SWING LATCHES
FOR EASY OPEN AND CLOSE

EXHAUST REQUIREMENTS:
220-340 CFM
STATIC PRESSURE:
.34-.67"

(4) HEAVY DUTY 8" CASTERS
(3) TOTAL LOCK
(1) DIRECTIONAL LOCK

HOODS SHOWN IN OPEN POSITION



HOODS SHOWN IN CLOSED POSITION

MOPEC LEGAL NOTE THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION OF MOPEC AND MAY NOT BE COPIED IN WHOLE OR PART, AND NO USE SHALL BE MADE OF IT OTHER THAN FOR WHICH IT HAS BEEN SUBMITTED WITHOUT THE WRITTEN PERMISSION OF A MOPEC MANAGER MOPEC 2026	UNLESS OTHERWISE SPECIFIED ALL DIMENSIONS IN INCHES TOLERANCES 0.00 ± 0.1 in 0.000 ± 0.05 in 0.0000 ± 0.001 in 0.00000 ± 0.0010 in ANGLES <15° BREAK ALL CORNERS 0.063 X 45° ALL FILLET RADIUS R0.063 (0.125 in) INSPECT OVAL DIMENSIONS	MATERIAL: STEEL	MOPEC GROUP Mortechnik MOPEC MORTTECH DRAWN TO 3RD ANGLE PROJECTION 	DRAWN BY DATE dollmayer 1/13/2026	PART NUMBER HBI410-B			
				CHECKED BY DATE 1/13/2026	DESCRIPTION DISSECTION CART, SINGLE VENTED W/ LID ASSIST			
						APPROVED BY DATE MMCLAIN 1/13/2026	SIZE A2 FILE NAME HBI410-B APPROVAL.dft	B
						SCALE 1:8	DO NOT SCALE IF IN DOUBT ASK	SHEET 1 OF 2



