

Single Body Carrier

Mopec's Single Body Carrier is designed for the secure, discreet transport of decedents. Constructed from durable materials with reinforced handles, it allows individuals to safely and efficiently carry a body while maintaining dignity and control during transfer.



JB026

SPECIFICATIONS

Construction

- Accommodates the JC026 body tray only (sold separately)
- (3) 5" heavy duty locking carrier wheels
- (1) 5" heavy duty directional lock wheel
- Curved saddles to hold trays in place during transportation

Materials

- 11ga 1.5" SQ 304 Stainless steel frame with a 180-grit polished finish, sheetmetal has a #4 finish
- (2) tray saddles to hold the tray

Dimensions

- 70.3" L X 25.3" W X 36" H

Capacity

- 375 lbs

Compatibility

- JC026 tray
- IR-series racks

Note:

- Trays sold separately

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI