

# Driven Hydraulic Forklift

Mopec's lifts are engineered to deliver effortless propulsion and smooth, controlled lifting. Designed for efficiency, safety, and maneuverability in mortuaries, laboratories, and medical facilities, these lifts enable operators to transport and elevate loads with precision and confidence.



JD405

## SPECIFICATIONS

### Construction

- 12V D.C. battery-powered hydraulic lift
- Built-in charger 12 V.D.C.
- (4) 5" drive wheels
- Electronic power drive system controlled by the throttle

### Capacity

- 1,000 lbs

### Materials

- Full steel construction and forks
- Hydraulic cylinder and hoses for controlling the lift
- Hard plastic controller and throttle

### Dimensions

- 44.5 L x 30.4" W x 78" H  
Fits through a standard 80" door
- Fork height: 12"-75" (min-max)
- 32" L forks (18" load center)

### Compatibility

- JD460 reversible forks (optional)  
Adds up to 19" of elevation

### Related Part Numbers

- JD450 64" hydraulic forklift with a max fork height of 64"
- Lift is 80" tall  
Fits through a standard 80" door
- JD451 76" hydraulic forklift with a max fork height of 76"
- Lift is 93" tall  
Does not fit through a standard 80" door

**MOPEC**  
GROUP

**Mopec**  
ELEVATING TECHNOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI