

Double Scrub Station, Optional ADA

Enhance hygiene, efficiency, and accessibility in your laboratory with our line of Scrub Stations, a durable and easy-to-maintain solution designed to support rigorous infection control standards. With an ADA-compliant model available, our Scrub Stations ensure inclusivity without compromising performance, making it an ideal choice for laboratories prioritizing safety, compliance, and user convenience.



SS102



SS102-ADA

SPECIFICATIONS

Construction

- Heliarc-welded seams and joints, ground and polished to match adjoining surfaces

Materials

- 304 stainless steel construction with a polished No. 4 finish

Utilities / Power

- 1/2" hot and 1/2" cold water supply lines
- 1 1/2" waste drain size
- (4) AA batteries required to power soap dispenser (OA-152)

Capacity

- SS102: 31 gallons of water
- SS102-ADA: 29.75 gallons of water

Performance

- Thermostatically-controlled mixing faucet with gooseneck spout
- Backflow protection for all plumbing fixtures to prevent backflow of contaminated water
- Rapid and positive drainage
- Access panels to get to plumbing
- ADA wheelchair accessible: (SS102-ADA)

Dimensions

- SS102: 63" L x 28" D x 36" H (48" H when station is mounted 12" above finished floor)
- SS102-ADA: 64" L x 23" D x 48" H
- SS102 sink depth: 14"
- SS102-ADA sink depth: 11"

Options & Accessories

- MB1046: Hot and cold foot valve (SS102 only)

Variations

| Model No. | Description | Dimensions | Sink Depth | Capacity |
|-----------|---------------------------|-----------------------|------------|---------------|
| SS102 | Double Scrub Station | 63" L x 28" D x 36" H | 14" Deep | 31 Gallons |
| SS102-ADA | Double Scrub Station, ADA | 64" L x 23" D x 48" H | 11" Deep | 29.75 Gallons |

MOPEC
GROUP

mopec
ELEVATING HYGIENE

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI