


# Triple Scrub Station, Optional ADA

Enhance hygiene, efficiency, and accessibility in your laboratory with our line of Scrub Stations, a durable and easy-to-maintain solution designed to support rigorous infection control standards. With an ADA-compliant model available, our Scrub Stations ensure inclusivity without compromising performance, making it an ideal choice for laboratories prioritizing safety, compliance, and user convenience.



SS103



SS103-ADA

## SPECIFICATIONS

### Construction

- Heliarc-welded seams and joints, ground and polished to match adjoining surfaces

### Materials

- 304 stainless steel construction with a polished No. 4 finish

### Utilities / Power

- 1/2" hot and 1/2" cold water supply lines
- 1 1/2" waste drain size
- (4) AA batteries required to power soap dispensers (OA-152)

### Capacity

- SS103: 51 gallons of water
- SS103-ADA: 31.7 gallons of water

### Performance

- Thermostatically-controlled mixing faucet with gooseneck spout
- Backflow protection for all plumbing fixtures to prevent backflow of contaminated water
- Rapid and positive drainage
- Access panels to get to plumbing
- ADA wheelchair accessible: (SS103-ADA)

### Dimensions

- SS103: 92" L x 28" D x 36" H (48" H when station is mounted 12" above finished floor)
- SS103-ADA: 96" L x 23" D x 48" H
- SS103 sink depth: 14"
- SS103-ADA sink depth: 11"

### Options & Accessories

- MB1046: Hot and cold foot valve (SS103 only)

## Variations

Model No.	Description	Dimensions	Sink Depth	Capacity
SS103	Triple Scrub Station	92" L x 28" D x 36" H	14" Deep	51 Gallons
SS103-ADA	Triple Scrub Station, ADA	96" L x 23" D x 48" H	11" Deep	31.7 Gallons

**MOPEC**  
GROUP

**mopec**  
ELEVATING HYGIENE

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI