





NOTES: UNLESS OTHERWISE SPECIFIED



1. APPLICABLE STANDARDS AND SPECIFICATIONS;
ASME Y14.5-2009, DIMENSIONING AND TOLERANCING
ASME Y14.38-2007, ABBREVIATIONS



CUSTOMER CUSTOMER NAME		APPROVAL		PROPRIETARY THIS DOCUMENT AND THE DATE DISCLOSED HEREIN AND HERewith IS NOT TO BE REPRODUCED, USED, OR OTHERWISE DISCLOSED IN WHOLE OR IN PART TO ANYONE WITHOUT WRITTEN AUTHORIZATION FROM MORTECH MFG.			 411 North Aerojet Avenue Azusa, CA 91702 tel: (800) 410-0100 (626) 334-1471 fax: (626) 334-1704 www.mortechmfg.com	
PROJECT PROJECT NAME		SIGNATURE		DATE	DO NOT SCALE DRAWING	APPROVALS DRAWN BY LUIS F. CHECKED BY APPROVED BY	DATE 01/16/19	TITLE MODEL 1036-240-72 VENTILATED HOOD DISSECTION TABLE
					MATERIAL 304 STAINLESS STEEL		SIZE C	DWG NO. 1036-240-72AP
					FINISH #4 SATIN		SCALE 1:12	REV. X01
							NEXT ASSY.	SHEET 1 OF 2

UTILITY REQUIREMENT NOTES:

1. FLUSH LINES CLEAR BEFORE CONNECTING TO STATION.
2. PLUMBER TO PROVIDE SHUT-OFF VALVES FOR WATER CONNECTION AND P-TRAP FOR DRAIN.
3. ELECTRICIAN TO PROVIDE LIQUID TIGHT JUNCTION BOX AND CONDUIT CONNECTIONS PER LOCAL CODES.
4. EXHAUST RECOMMENDATIONS: PROVIDE 450-400 cfm @ 0.50w.g.



<p>CUSTOMER CUSTOMER NAME</p> <p>PROJECT LINE 1 LINE 2</p>	<p>APPROVAL</p> <p>SIGNATURE _____ DATE _____</p>	<p style="text-align: center;">PROPRIETARY</p> <p style="text-align: center; font-size: small;">THIS DOCUMENT AND THE DATE DISCLOSED HEREIN AND HEREWITH IS NOT TO BE REPRODUCED, USED, OR OTHERWISE DISCLOSED IN WHOLE OR IN PART TO ANYONE WITHOUT WRITTEN AUTHORIZATION FROM MORTECH MFG.</p>	<p style="text-align: right; font-size: small;">MOPEC GROUP 411 North Aerojet Avenue Azusa, CA 91702 tel: (800) 410-0100 (626) 334-1471 fax: (626) 334-1704 www.mortechmfg.com</p> <p style="text-align: center;">MODEL 1036-240-72 ROUGH IN</p>
		DO NOT SCALE DRAWING	TITLES
		MATERIAL 304 STAINLESS STEEL	APPROVALS
		FINISH #4 SATIN	DATE
		DRAWN BY LUIS F.	DWG NO. 1036-240-72AP
		CHECKED BY	REV. X01
		APPROVED BY	SCALE 1:8
		NEXT ASSY.	SHEET 2 OF 2