




NOTES: UNLESS OTHERWISE SPECIFIED


1. APPLICABLE STANDARDS AND SPECIFICATIONS;
ASME Y14.5-2009, DIMENSIONING AND TOLERANCING
ASME Y14.38-2007, ABBREVIATIONS



CUSTOMER CUSTOMER NAME		APPROVAL		PROPRIETARY THIS DOCUMENT AND THE DATE DISCLOSED HEREIN AND HEREWITH IS NOT TO BE REPRODUCED, USED, OR OTHERWISE DISCLOSED IN WHOLE OR IN PART TO ANYONE WITHOUT WRITTEN AUTHORIZATION FROM MORTECH MFG.			MOPEC GROUP 411 North Aerojet Avenue Azusa, CA 91702 tel: (800) 410-0100 (626) 334-1471 fax: (626) 334-1704 www.motechmfg.com	
PROJECT PROJECT NAME		SIGNATURE		DATE	DO NOT SCALE DRAWING	APPROVALS	DATE	TITLE MODEL 1036-240-84
					MATERIAL 304 STAINLESS STEEL	DRAWN BY LUIS F.	01/16/19	VENTILATED HOOD DISSECTION TABLE
					FINISH #4 SATIN	CHECKED BY		SIZE C
						APPROVED BY		DWG NO. 1036-240-84AP
								REV. X01
								SCALE 1:12
								NEXT ASSY.
								SHEET 1
								OF 2

UTILITY REQUIREMENT NOTES:

1. FLUSH LINES CLEAR BEFORE CONNECTING TO STATION.
2. PLUMBER TO PROVIDE SHUT-OFF VALVES FOR WATER CONNECTION AND P-TRAP FOR DRAIN.
3. ELECTRICIAN TO PROVIDE LIQUID TIGHT JUNCTION BOX AND CONDUIT CONNECTIONS PER LOCAL CODES.
4. EXHAUST RECOMMENDATIONS: PROVIDE 450-400 cfm @ 0.50w.g.



CUSTOMER CUSTOMER NAME		APPROVAL		PROPRIETARY			 <small>411 North Aerojet Avenue Azusa, CA 91702 tel: (800) 410-0100 (626) 334-1471 fax: (626) 334-1704 www.mortechmfg.com</small>	
PROJECT LINE 1 LINE 2		SIGNATURE		<small>THIS DOCUMENT AND THE DATE DISCLOSED HEREIN AND HEREWITH IS NOT TO BE REPRODUCED, USED, OR OTHERWISE DISCLOSED IN WHOLE OR IN PART TO ANYONE WITHOUT WRITTEN AUTHORIZATION FROM MORTECH MFG.</small>			MODEL 1036-240-84 ROUGH IN	
		DATE		DO NOT SCALE DRAWING			TITLE	
				MATERIAL 304 STAINLESS STEEL			1036-240-84AP	
				FINISH #4 SATIN			SIZE C DWG NO. 1036-240-84AP REV. X01	
				APPROVED BY			SCALE 1:8 NEXT ASSY. SHEET 2 OF 2	