

Autopsy Carrier Chassis, Elevating

Mopec carriers are designed to provide safe, secure, and dignified transport in professional environments. Built from durable, high-quality materials, they offer dependable strength and stability while supporting efficient handling. Thoughtful design features help improve ease of use, making Mopec carriers a reliable choice for everyday professional needs.



600018-FRC-JC023



600018-FRC-JC027

SPECIFICATIONS

Construction

- Fits corresponding JC Series 23" or 27" trays (sold separately)
- Front and rear hand cranking mechanism to control up to 6" elevation and tilt
- (4) 8" Casters with brakes
- Cadaver lift accessible

Materials

- 11 Gauge 304 stainless steel

Dimensions

- 74.9" L x 26" W x 30.4" (minimum) – 36.4" (maximum) H

Capacity

- 600018-FRC-JC023: 300 lbs
- 600018-FRC-JC027: 350 lbs

Performance

- Hand crank one side independently, or both at same time

Compatibility

- 600018-FRC-JC023: JC023 and JC023D trays only
- 600018-FRC-JC027: JC027 and JC027D trays only

Note:

- Trays sold separately.

Variations

Model No.	Description	Dimensions
600018-FRC-JC023	Elevating Carrier Chassis for JC Series 23" Trays	74.9" L x 26" W x 30.4" - 36.4" H
600018-FRC-JC027	Elevating Carrier Chassis for JC Series 27" Trays	74.9" L x 26" W x 30.4" - 36.4" H

MOPEC
GROUP



www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI