





NOTES: UNLESS OTHERWISE SPECIFIED

1. APPLICABLE STANDARDS AND SPECIFICATIONS;
ASME Y14.5-2009, DIMENSIONING AND TOLERANCING
ASME Y14.5-2007, ABBREVIATIONS

TRAY CHANNEL WITH DRAIN HOLE
304 SS, 3" DIA DRAIN

CARRIER ACCEPTS MOPEC TRAY
MODELS: JC027
UNDERSIDE OF TRAY MATES TO THE TOP FRAME
OF CARRIER TO SECURELY HOLD TRAY IN PLACE
WITHOUT MECHANICAL FASTENERS

DESIGN LOAD: 500 LBS.



CUSTOMER		APPROVAL		PROPRIETARY		MOPEC GROUP		411 North Aerojet Avenue Azusa, CA 91702 tel: (800) 410-0100 (626) 334-1471 fax: (626) 334-1704 www.mortechmfg.com	
PROJECT				THIS DOCUMENT AND THE DATE DISCLOSED HEREIN AND HERewith IS NOT TO BE REPRODUCED, USED, OR OTHERWISE DISCLOSED IN WHOLE OR IN PART TO ANYONE WITHOUT WRITTEN AUTHORIZATION FROM MORTECH MFG.		TITLE		MODEL 600018-JC027	
LINE 1 LINE 2		SIGNATURE		DATE		DO NOT SCALE DRAWING		APPROVALS	
						MATERIAL		DRAWN BY	
						FINISH		CHECKED BY	
						GRAINED		APPROVED BY	
						SIZE		DWG NO.	
						C		600018-JC027	
						SCALE		NEXT ASSY.	
						1:8		SHEET	
								1 OF 1	
								REV. X01	