

Autopsy Tray Carrier

Mopec carriers are designed to provide safe, secure, and dignified transport in professional environments. Built from durable, high-quality materials, they offer dependable strength and stability while supporting efficient handling. Thoughtful design features help improve ease of use, making Mopec carriers a reliable choice for everyday professional needs.



600018-JC023



600018-JC027

SPECIFICATIONS

Construction

- (4) 8" heavy duty casters with brake
- Tray channels hold the tray in place during cart movement
- 1.5" slope to aid in drainage

Materials

- 11-Gauge 304 stainless steel

Dimensions

- 600018-JC023: 74.9" L X 19.9" W X 32.8" to 34.2" H
Inner tray channel width: 18.8"
- 600018-JC027: 74.9" L X 23.8" W X 32.8" to 34.2" H
Inner tray channel width: 22.8"

Capacity

- 600018-JC023: 350 lbs
- 600018-JC027: 375 lbs

Compatibility

- 600018-JC023: JC023 and JC023D trays only
- 600018-JC027: JC027 and JC027D trays only

Note:

- Trays sold separately.

Variations

Model No.	Description	Dimensions
600018-JC023	Autopsy Carrier for JC Series 23" Trays	74.9"L x 19.9"W x 32.8"-34.2"H
600018-JC027	Autopsy Carrier for JC Series 23" Trays	74.9"L x 23.8"W x 32.8"-34.2"H

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI