


# Autopsy Tray

Mopec's body and autopsy trays are engineered from high-quality stainless steel to withstand daily use in pathology and mortuary settings. Each tray reflects Mopec's commitment to durable, professional-grade equipment that meets the demands of clinical and post-mortem environments.



GA031

## SPECIFICATIONS

### Construction

- Features a flat surface with raised edges
- Long sides of the tray are sloped to allow for easy loading and unloading
- Creased to allow for draining to the included drain hole, plug and chain

### Materials

- 304 stainless steel, 14GA sheetmetal has a polished No. 4 finish

### Dimensions

- 77.5" L x 31" W x 1.5" H

### Compatibility

- IA32/IB32/IJ32 Racks
- JD5001 Scissor Lift
- DA031 Autopsy Cart Chassis

### Capacity

- 500 lbs

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI