


I 2 3 4 5 6 7 8 9 10 II

REV	DESCRIPTION	DATE	ENR #	ISS
A	INITIAL RELEASE	00-XXXX-19	XXXXX	MM



Item #	Qty	File Name	Title	Material	Rev
1	1	HO1112-B.asm	ELEVATING CART FRAME ASSEMBLY		B
2	8	PB0016.par	NUT - SS, 1/4"-20, NYLON LOCK		A
3	8	PB0137.par	BOLT - ZINC, M8 x 1.25 x 14mm, HEX HEAD		A
4	8	PB0211.par	WASHER - SS, M8 x 5/16", SERRATED BELLEVILLE		A
5	8	PB0303.par	BOLT - SS - 1/4"-20 x 2-1/4" PARTIAL THREAD		A
6	16	PB0603.par	BOLT HEX M5 X 45MM SS FULL THREAD	McMASTER-CARR # 91287A331	A
7	3	PM0803.par	5"X1" CASTER WHEEL, WITH BRAKE, GREY LOCK		A
8	1	PM0804.par	5"X1" CASTER WHEEL, WITH BRAKE, DRIVE LOCK, GREEN		A
9	1	PY0214.asm	LIFT - SUSPA 4 CYL 1000 LBS 12" STROKE HAND CRANK		A
10	1	SB141-A.asm	ELEVATING CART TOP ASSEMBLY		A

**MOPEC LEGAL NOTE**  
 THIS DOCUMENT CONTAINS  
 PROPRIETARY INFORMATION OF  
 MOPEC, AND MAY NOT BE COPIED IN WHOLE  
 OR PART, AND NO USE SHALL BE MADE OF  
 IT OTHER THAN FOR WHICH IT HAS BEEN  
 SUBMITTED WITHOUT THE WRITTEN  
 PERMISSION OF A MOPEC MANAGER

UNLESS OTHERWISE SPECIFIED  
 ALL DIMENSIONS IN INCHES  
 TOLERANCES 0.0 ± 0.1 in  
 0.00 ± 0.05 in  
 0.0000 ± 0.0010 in  
 0.0000 ± 0.0010 in  
 ANGLES ± 0.5°  
 BREAK ALL CORNERS 0.063 X 45°  
 ALL FILLET RADIUS R0.063 (0.125 in)  
 INSPECT OVAL DIMENSIONS

MATERIAL:  
 UNLESS OTHERWISE SPECIFIED  
 ALL DIMENSIONS IN INCHES

**MOPEC GROUP**  
 Mortechnik  
 MOPEC GROUP

DRAWN BY: dalmayer  
 DATE: 1/23/2026  
 CHECKED BY:  
 DATE: 1/23/2026  
 APPROVED BY:  
 DATE: MMCCCLAIN 1/23/2026


PART NUMBER	HA230-A
DESCRIPTION	CART, DISSECTION, HYDRAULIC HAND CRANK
SIZE	A2
FILE NAME	HA230-A.dft
SCALE	1:5.56
DO NOT SCALE IF IN DOUBT, ASK	
SHEET 1 OF 2	

H

H

I 2 3 4 5 6 7 8 9 10 II

\*HA230-A.dft\*



REV	DESCRIPTION	DATE	ENR #	ISS
A	INITIAL RELEASE	00-XXXX-19	XXXXXX	MM

**MOPEC LEGAL NOTE**  
 THIS DOCUMENT CONTAINS  
 PROPRIETARY INFORMATION OF  
 MOPEC, AND MAY NOT BE COPIED IN WHOLE  
 OR PART, AND NO USE SHALL BE MADE OF  
 IT OTHER THAN FOR WHICH IT HAS BEEN  
 SUBMITTED WITHOUT THE WRITTEN  
 PERMISSION OF A MOPEC MANAGER  
 MOPEC 2026

UNLESS OTHERWISE SPECIFIED  
 ALL DIMENSIONS IN INCHES  
 TOLERANCES 0.0 ± 0.1 in  
 0.00 ± 0.05 in  
 0.000 ± 0.010 in  
 0.0000 ± 0.0010 in  
 ANGLES ± 0.5°  
 BREAK ALL CORNERS 0.063 x 45°  
 ALL FILLET RADI R0.063 (0.125 in)  
 INSPECT OVAL DIMENSIONS

MATERIAL:  
 MOPEC GROUP  
 MOPEC MORTTECH  
 DRAWN TO SRED ANGLE PROJECTION

MOPEC GROUP  
 MOPEC MORTTECH  
 DRAWN BY: dalmayer  
 DATE: 1/23/2026  
 CHECKED BY:  
 DATE: 1/23/2026  
 APPROVED BY:  
 DATE: MMCCCLAIN 1/23/2026

PART NUMBER: HA230-A  
 DESCRIPTION: CART, DISSECTION, HYDRAULIC HAND CRANK  
 SIZE: A2  
 FILE NAME: HA230-A.dft  
 SCALE: 1:6.67  
 DO NOT SCALE IF IN DOUBT ASK  
 SHEET 2 OF 2

\*HA230-A.dft\*