

# Cadaver Lift, End Access

Mopec's lifts are engineered to deliver effortless, smooth, and controlled lifting. Designed for safe maneuverability in mortuaries, laboratories, and medical facilities. These lifts enable operators to transport and elevate loads with precision and confidence.



JD5001

## SPECIFICATIONS

### Construction

- End loading only
- (2) full side rails and two end stop gates
- (4) 5" Full-swivel casters with brakes
- Integrated lift controller with scale display

### Capacity

- 500 lbs (tray weight limit)

### Materials

- 11GA 304 Stainless steel tubing, various gauges of stainless steel sheet metal
- Galvanized steel rollers
- Stainless steel transfer balls

### Dimensions

- 95.2" L x 34.4" W (36.4" including side handle) x 13.4" (minimum) – 76.8" (maximum) H

### Utilities / Power

- Battery-operated electromechanical, scissor lifting system
- 24V power system, (2 Rechargeable deep-cycle 12-volt batteries included) battery cord accepts any standard 120V/1hp/60V outlet

### Compatibility

- JC Series 23" and 27" trays
- GA100 32" Legacy tray
- Note: Only compatible with Mopec trays, not compatible with Mortech Legacy trays

### Note:

- Trays sold separately.

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI