

# Ventilated Hood Dissection Table

Mopec's Ventilated Hood Dissection Tables are designed to provide a durable and functional workspace for anatomy labs and pathology settings. Built for reliable daily use, it offers a stable work surface and practical design to support efficient specimen handling and examination. Ideal for educational and clinical environments, this table is constructed to meet the demands of professional laboratory use.



1036-240-60



1036-240-84



1036-240-72

## SPECIFICATIONS

### Construction

- TÜV Certified
- Heliarc-welded seams and joints, ground and polished to match adjoining surfaces
- Radii on all inside and outside corners for easy cleaning
- Stainless steel base frame with adjustable feet
- Faucet with cold water control; hose rinse assembly with cold water control
- Hinged plexiglass work shield
- (1) LED light fixture and (1) flex-arm LED light
- GFCI duplex receptacle with waterproof cover
- (2) Access panels (one per side)
- Plumbing and electrical lines are factory-installed and require only single-point connection at the project location

### Utilities / Power

- Water: Provide 60–65 PSI water pressure and 4 GPM of flow
- 1/2" Hot, 3/4" Cold water and 1 1/2" drain rough-ins
- Provide 450–400 CFM at 0.50 W.G 6" flex duct connection
- Electrical: Provide (1) circuits in 4" x 4" junction box: 115V/60Hz / 1Ph / 15A

### Materials

- 304 Stainless steel with No. 4 finish

### Dimensions

- 30" L X 56" H
- 1036-240-60: 60" W
- 1036-240-72: 72" W
- 1036-240-84: 84" W
- Sink: 12" L X 14" W X 8" D

### Variations

Model No.	Description	Dimensions
1036-240-60	Ventilated Hood Dissection Table	30"L x 60"W x 56"H
1036-240-72	Ventilated Hood Dissection Table	30"L x 72"W x 56"H
1036-240-84	Ventilated Hood Dissection Table	30"L x 84"W x 56"H

**MOPEC**  
GROUP

**mopec**  
CLEANING TECHNOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI