

# Wall Mount Service Station

Mopec's Water Control Units are a durable, stainless steel water management panel designed for efficient hot and cold water regulation in laboratory and pathology workstations. Engineered for reliability and ease of use, these models feature a compact design for seamless integration with lab equipment. Built to withstand demanding environments, it supports safe, sanitary, and streamlined workflow operations.



1036-2



1036-2SE

## SPECIFICATIONS

### Construction

- TÜV Certified
- Fully fabricated stainless steel
- Spray hose assembly with cold water control valve, nozzle and 10 ft./3 m flexible hose
- Waste disposal unit
- Hot and cold-water tempering unit with lever handle
- Utility faucet with elbow handle and swivel spout
- GFCI duplex receptacle
- Deluxe aspirator with reverse flow for additional waste removal: "reverse flow" feature instantly relieves suction clogging by providing a reverse pressure of water
- Open stainless-steel shelving designed for storage of chemical bottles
- All plumbing and electrical hook-ups designed for easy installation

### Utilities / Power

- Water – Provide 60–65 PSI water pressure and 4 GPM of flow.
- Electrical – Provide 4" x 4" junction box 115V / 60 Hz / 1 Ph / 20A

### Materials

- 304 Stainless steel with #4 finish

### Dimensions

- 30" L X 24.5" W X 71" H
- 1036-2SE: 85.75" H with shower

### Options & Accessories

- 1036-2SE includes emergency shower and eyewash assemblies only.

## Variations

Model No.	Description	Dimensions
1036-2	Wall Mount Service Station	30"L x 24.5"W x 71"H
1036-2SE	Wall Mount Service Station with Shower/Eyewash	30"L x 24.5"W x 85.75"H

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI