

# Mobile Scale Stand

Mopec's Mobile Scale Stand is a durable stainless steel wheeled stand designed to support compatible autopsy and pathology scales. It features an offset center post and a weighted base to provide stability during weighing procedures. Four heavy-duty locking casters allow the stand to be easily moved and securely positioned within the workspace.



BB035

## SPECIFICATIONS

### Construction

- Mobile stand
- The offset center post and weighted base provide secure operation
- Durable and easy to clean
- (4) industrial locking casters provide complete mobility

### Materials

- Type 304 stainless steel with #4 finish

### Dimensions

- 22.5" L X 24" W (Base) X 84.0" H

### Compatibility

- BB002: Autopsy Scale with Bow, Ring and Perforated Pan
- BB003: Autopsy Scale with Bow, Ring and Pan
- BB013: Autopsy Scale with Bowl
- BB014: Autopsy Scale with Bowl with Drain Holes

### Options & Accessories

- BB007: Autopsy Scale Perforated Pan
- BB009: Autopsy Scale Solid Pan
- BB037: Autopsy Scale Bow and Ring Set with Perforated Pan

### Note:

- **Compatibility and Options & Accessories products are sold separately.**

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI