

# Down-Back Draft Grossing Station

Mopec's Down-Back Draft Grossing Station is designed to support a safer and more comfortable environment during surgical grossing procedures. The workstation helps manage fumes and odors at the work area while maintaining an ergonomic, user-focused design. Built for durability and everyday use, it is well suited for the demands of busy pathology labs.



BF550

## SPECIFICATIONS

### Construction

- Facility ventilation pulls between 300-350CFM
- 60" High pressure hand spray hose and 1.5" drainpipe
- Top mounted, watertight G.F.C.I. Outlet
- 18" Standard Fluorescent Replacement Bulb with Lighted On/Off Switch

### Materials

- 18 Gauge 304 Stainless Steel with a No. 4 Finish casing
- 16 Gauge 304 Stainless Steel with a No.4 Finish frame

### Dimensions

- 21" L X 36" W X 58.7" H
- Workstation is at 37.3" H
- 18" x 4" rectangular duct

### Utilities / Power

- Hospital Grade Line Cord 115v/1ph/60Hz (8' long)

### Options & Accessories

- BF500 (recirculating version of this table)
- BF003 Particle Filter (1 included)
- BF005 Potassium Permanganate Filter (1 included)
- BF101 Replacement Pad (1 pack of 3 included)

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI