

# Formalin Dispensing Cube

This formalin dispensing cube provides a safe and efficient solution for storing and dispensing formalin in laboratory and pathology environments. Designed for controlled handling, it helps minimize exposure to fumes while allowing for precise, convenient dispensing. Built for durability and ease of use, it supports a cleaner, safer workspace and consistent workflow.



FP025

## SPECIFICATIONS

### Construction

- Waterproof on/off switch
- Formalin dispensing spigot with adjustable flow rate
- 6' C13 Hospital-Grade Power Cord
- Pump: 1.5 gpm, built in 45psi internal bypass, 0.54amps at 60 psi, 115vac 50 -100 psi, 6GPM
- Connects to standard formalin cubitainer
- Screw cap on your cubitainer, then click on quick-disconnect fitting to dispense

### Dimensions

- 17.1" L X 8.4" W X 15.4" H

### Utilities / Power

- Temperature Range: 70.0±10.0°F or 21.0±6.0°C
- Relative Humidity Range: RH 35%-75%
- Elevation Range: 0-2000m
- Electrical Specification: 115v/1ph/60hz

### Materials

- 14/18 Gauge 304 stainless steel with #4 finish

### Options & Accessories

- BE092: Replacement Pad
- BG520: Mopec Prefilled Formalin Cube, 2.5 Gallons
- BG525: Mopec Prefilled Formalin Cube, 5 Gallons

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI