

Bucket Holder, ADA

Mopec has various dissection accessories that are built for seamless integration, they help optimize workspace utilization, increase safety, and provide added flexibility for different applications and user preferences.



H0009

SPECIFICATIONS

Construction

- Designed to hold BA005 bucket on the exterior of the cart
- Will not interfere with hoods or ADA requirements
- Includes a locking pin

Materials

- 304 stainless steel construction
- 1" SQ. x 11 GA tubing

Dimensions

- 2.5" (arm-in) - 9.2" W (arm out) x 9.1" H

Performance

- Is removable (bolt-on)
- Can swivel when not in use

Compatibility

- BA005
- HA230-ADA
- HB1550

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI