

Large Animal Trimming Table with Casters

Mopec's Large Animal Trimming Table features a durable stainless-steel, kidney-shaped work surface designed to improve operator access and specimen handling during large animal trimming procedures. Its spacious design supports efficient workflow, while heavy-duty locking casters provide smooth mobility and stable positioning within veterinary, pathology, and research lab environments.



LA100

SPECIFICATIONS

Construction

- (4) Extreme-duty 6" Caster locking casters
- Scissor lift for raising and lowering the large kidney shaped table
- Equipped with accordion cover for wash-down to protect lift system
- Handheld elevating switch for raising and lowering the scissor lift
- 8' long 115V/ 1hp/ 60hz cord for powering the scissor lift

Materials

- 11 GA 304 stainless steel top with a No.4 finish
- Electro-Hydraulic powered scissor lift
- 0.25" thick 304 stainless steel tube base

Dimensions

- 84" L X 48" W X 17.5" to 53.5" H

Capacity

- 4000 lbs
- 2000 lbs if offset

Utilities / Power

- 115V / 1hp / 60hz electrical requirement for powering the scissor lift

MOPEC
GROUP

Mopec
ELEVATING TECHNOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI