

# Stainless Steel Scullery Sink, 48"

Mopec's Scullery Sinks are a heavy-duty stainless steel sink designed for demanding mortuary and laboratory environments. Its deep basin and durable construction support thorough cleaning and decontamination tasks, while the seamless design promotes easy sanitation. Built for reliability and efficiency, it provides a practical solution for maintaining a clean, professional workspace.



MS2448LR-24



MS2448L-24



MS2448R-24

## SPECIFICATIONS

### Construction

- Hot and cold water fixtures with gooseneck faucet and wrist blade handles
- Rounded sink and drainboard
- Leg diameter: 1.6"
- Drain diameter: 3.6"
- Drainboard to left or right of sink, or on both sides
- (4) Leveling feet

### Materials

- 304 stainless steel, 16 GA

### Dimensions

- MS2448L-24 and MS2448R-24: 74" L X 29.4" W X 45" H
- MS2430LR-24: 96" L X 29.4" W X 45" H
- Sink: 48" L X 24" W X 14" H

### Utilities / Power

- Provide 60-65 PSI water pressure and 4 GPM of flow

### Options & Accessories

- Left, right, or double-sided options available.

## Variations

Model No.	Description	Dimensions
MS2430L-24	Scullery Sink, Drainboard on Left	74"L x 29.4"W x 45"H
MS2430R-24	Scullery Sink, Drainboard on Right	74"L x 29.4"W x 45"H
MS2430LR-24	Scullery Sink, Drainboard on Both Sides	96"L x 29.4"W x 45"H



www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI