

# Vacuum Valve Deck Mount

Mopec's Deck-Mounted Vacuum Valve provides a reliable solution for controlling vacuum lines in laboratory and clinical environments. Designed for direct installation onto work surfaces, it offers convenient access while supporting an organized and efficient workspace. Built for durability and consistent performance, it ensures safe and dependable vacuum management in demanding settings.



OA-102-DM

## SPECIFICATIONS

### Construction

- Deck-Mounted Laboratory vacuum valve with turn on/off valve and serrated hose fitting.
- Mounted to the deck of the workstation
- Max pressure: 75 PSI
- Lever activated air flow

### Dimensions

- 5.8" L X 3.9"H X 1.9" (Base Diameter)
- 3/8" inlet

### Materials

- Male 3/8" forged brass serrated inline nozzle with lever activation

### Environmental Conditions

- Max Temperature: 52°C
- Min. Temperature: 0°C

**MOPEC**  
GROUP

**Mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI