

Eyewash Assembly, Hose Style

Mopec's Hose-Style Eyewash Assembly provides a flexible and effective solution for emergency eye rinsing in laboratory and clinical environments. Designed for quick response, the handheld hose allows for targeted flushing, making it easy to direct water where it's needed most. Built for reliability and ease of use, it supports workplace safety.



OA-111

SPECIFICATIONS

Construction

- Hand-held
- 3.2 gpm once activated
- Single spray head
- Squeeze handle to begin the flow of water
- Release handle to stop the flow of water

Dimensions

- 1 – 1/2" Spray head diameter
- 3/8" inlet pipe size

Materials

- 8' PVC hose
- HDPE spray head

Utilities / Power

- 300 PSI maximum water pressure

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI