

Eyewash Wall Mount with 90° Swivel

Mopec's OA-113 wall-mounted eyewash features a durable, space-saving design with a 90-degree swivel for flexible positioning and easy access in emergency situations. Engineered for reliable performance, it enables rapid eye irrigation to help flush contaminants and reduce injury risk in laboratory and clinical environments.



OA-113

SPECIFICATIONS

Construction

- 1.4 gpm once activated
- Activates water flow when pulled down 90°
- Stops water flow when lifted 90°
- (2) Hand spray heads so each eye has dedicated eyewash

Materials

- Stainless steel
- HDPE spray heads

Dimensions

- 14.9" L X 8.2" W X 4" H
- 0.5" faucet tubing

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI