

# Total Draft Trimming Table with Sink

Mopec's Total Draft Trimming Table is a specimen trimming workstation with integrated downdraft ventilation to help remove fumes and improve air quality. It features a built-in sink, perforated work surface for drainage, and electrical outlets, providing a durable and hygienic solution for pathology and laboratory environments.



## SPECIFICATIONS

### Construction

- Grid plate work surfaces with 1/2" recessed holes to allow for drainage
- 10' Spray hose assembly that can comfortably reach the length of the table
- Adjustable airflow at in-house ventilation hook-up

### Materials

- 14 GA 304 Stainless steel with No. 4 finish
- Pex piping for plumbing purposes

### Dimensions

- 96" L X 30" W X 37" H
- Downdraft Bellows: 60.5" H
- Small Sink: 11" L X 20" W X 10" D
- Large Sink: 62" L X 20.5" W

### Capacity

- 300 lbs

### Utilities / Power

- (2) G.C.F.I 120V/ 1hp/ 60hZ outlets with waterproof covers
- 1/2" Hot water, 1/2" Cold Water, and 1-1/2" Drain for water hookup
- House Ventilation capable of 570CFM at 0.70" WG

**MOPEC**  
GROUP

**mopec**  
ELEVATING PATHOLOGY

**MORTECH**  
MANUFACTURING

[www.mopecgroup.com](http://www.mopecgroup.com)

800.362.8491

800 Tech Row, Madison Heights, MI