


Total Draft Elevating Trimming Table with Sink

Mopec's Total Draft Elevating Trimming Table is a specimen trimming workstation with integrated downdraft ventilation to help remove fumes and improve air quality. It features a built-in sink, perforated work surface for drainage, and electrical outlets, providing a durable and hygienic solution for pathology and laboratory environments.



SPECIFICATIONS

Construction

- Grid plate work surfaces with 1/2" recessed holes to allow for drainage
- 10' Spray hose assembly that can comfortably reach the length of the table
- Adjustable airflow at in-house ventilation hook-up
- Adjust the height of the sink using an Up/Down switch

Materials

- 14 GA 304 Stainless steel with No. 4 finish
- Pex piping for plumbing purposes
- (2) 12" linear actuators to raise/lower the sink

Dimensions

- 78" L X 30" W X 32" to 44" H
- Duct stubs: 57.4" H

Capacity

- 300 lbs

Utilities / Power

- (2) G.C.F.I 120V/ 1hp/ 60hZ outlets with waterproof covers
- 1/2" Hot water, 1/2" Cold Water, and 1-1/2" Drain for water hookup
- House Ventilation capable of 560CFM at 0.60" WG

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI