


Hydraulic Autopsy Carrier

Mopec's Autopsy Carrier with Hydraulic Lift combines lightweight, corrosion-resistant construction with powerful, smooth hydraulic height adjustment. Its design is easy to clean and durable, while the hydraulic system enables effortless positioning for safe, efficient body transport. Designed for reliability and ease of use, it supports consistent performance in mortuary and pathology settings.



600017-H

SPECIFICATIONS

Construction

- Hydraulic Stryker base able elevate tray 14" and tilt tray according to preference
- Enable fifth wheel steering, neutral, or full brakes with foot controls
- For use with T3603 Plastic Body tray only (not included)
- Underside of T3603 tray mates to top of the carrier to securely hold the tray in place without the use of mechanical fasteners

Materials

- 11ga 1.5" SQ 304 Stainless steel frame with a 180-grit polish finish

Dimensions

- 80.5" L X 30" W X 26.8" to 40.1" H (lower to higher end)

Capacity

- 500 lbs

Performance

- Raise or lower the cart using the dual sided mounted foot controls

Compatibility

- T3603 (Plastic Body Tray)

Note:

- Trays sold separately.

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI