


Hydraulic Covered Cadaver Carrier with Rollers

Mopec's Hydraulic Autopsy Carrier with Rollers combines lightweight, corrosion-resistant construction with powerful, smooth hydraulic height adjustment. Its design is easy to clean and durable, while the hydraulic system enables effortless positioning for safe, efficient body transport. Designed for reliability and ease of use, it supports consistent performance in mortuary and pathology settings.



600040-H

SPECIFICATIONS

Construction

- Fits 23" JC series trays (sold separately)
- Raise or lower the cart 13" using the dual-sided mounted foot controls
- Enable fifth wheel steering, neutral, or full brakes with foot controls
- (4) 8" diameter wheels
- Removable aluminum tubing canopy and black vinyl cover
- (5) galvanized steel rollers for smooth transfer of the tray
- (2) locking latch mechanisms at the head and foot of the carrier

Materials

- 11 Gauge square 304 stainless steel carrier frame with a No. 4 polished finish
- Galvanized steel rollers
- The false frame is made of lightweight aluminum
- Hydraulic Stryker base
- Vinyl cover

Capacity

- 500 lbs

Dimensions

- 81" L X 26.5" W X 22" to 35" (50.7" to 63.7" H with false frame)
- Cover: 82" L X 33" W X 21" D

Compatibility

- JC023 and JC023D trays (included and available individually)
- JZ029: Black vinyl cover (included and available individually)

Note:

- Trays sold separately

MOPEC
GROUP

mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI