

Dissecting Cart, One-Sided Tilting

Mopec's Dissection Carts are durable, stainless steel workstations designed to support efficient and hygienic specimen handling in laboratories and educational settings. Built for mobility and stability, they feature a spacious work surface, integrated drainage, and smooth-rolling casters for easy maneuverability.



600042-A

SPECIFICATIONS

Construction

- Height-adjustable top can be tilted 6" from 35" to 41" off the floor
- (4) 8" locking casters
- 2" wide rim around the edges with a 1.5" deep recessed trough top
- 7.5" wide drainage spout with notch for hanging a bucket, or drain into an autopsy sink
- 19.5" square shelf

Materials

- 14-gauge 304 Stainless steel top with a #4 finish
- 11-gauge 304 Stainless steel frame with 180 grit finish
- 16 gage 304 Stainless steel shelf with #4 finish

Dimensions

- 35" L X 23" W X 35"- 41" H (adjustable height top)

Capacity

- 140 lbs

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI