

Dissecting Cart, Stationary Height

Mopec's Dissection Carts are durable, stainless steel workstations designed to support efficient and hygienic specimen handling in laboratories and educational settings. Built for mobility and stability, they feature a spacious work surface, integrated drainage, and smooth-rolling casters for easy maneuverability.



600042

SPECIFICATIONS

Construction

- Stationary top
- (4) 8" locking casters
- 2" wide rim around the edges with a 1.5" deep recessed trough top
- 7.5" wide drainage spout with notch for hanging a bucket, or drain into an autopsy sink
- 19.5" square shelf

Materials

- 14-gauge 304 Stainless steel top with a #4 finish
- 11-gauge 304 Stainless steel frame with 180 grit finish
- 16 gage 304 Stainless steel shelf with #4 finish

Dimensions

- 35" L X 23" W X 35" H

Capacity

- 160 lbs

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI