


100 mL Urine Collection and Transport Bottle

These 100 ml Urine Specimen Collection and Transport Bottles are pre-sealed with an internal cap integrity seal and feature clear volume markings. These bottles are available with or without urine temperature monitoring thermometers for accurate and secure specimen handling.



FS003 and FS004

SPECIFICATIONS

Construction

- Pre-sealed with an internal bottle cap integrity seal
- Volume markings
- Available with or without a urine temperature monitoring thermometer affixed to the outside of the cup

Capacity

- 100 ml

Quantity

- 250 per case

Performance

- Transport and collect urine samples

Variations

Model No.	Description
FS003	Without Strip
FS004	With Strip

MOPEC
GROUP

Mopec
ELEVATING PATHOLOGY

MORTECH
MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI